

THE INTELLIGENT CHOICE IN FALL PREVENTION



Smart new features and our Lifetime Warranty put Posey's new Sitter II™ fall monitor at the top of the class.

The new Sitter II restraint-free fall monitor incorporates the best features of the original Sitter monitor and adds exciting new capabilities that make it the intelligent choice for your fall prevention needs. Sitter II uses ten different alarm tones, including a selection of easily recognizable musical tunes, to alert caregivers when a patient attempts to get up from bed, a chair, wheelchair, or geri-chair.

The Sitter II offers a voice option which allows a family member, friend or caregiver to record a personalized message cautioning the patient against getting up unattended. With Sitter II, you can choose the voice option alone, the voice option with an audible tone, a tone only, or mute (using the nurse call function).

To save your staff time and trouble, the Sitter II offers a hold feature which allows the patient to be away from the fall monitor for up to 6 hours without the need to turn off the alarm and then reset it. The alarm is automatically reactivated when the patient returns to bed or chair.

The new Sitter II fall monitor is backed by Posey's unique Lifetime Warranty, making it a smart choice for your facility.

Sitter II
A smart choice in
fall prevention



POSEY *CA*
Care Alternatives Division

The Posey® Sitter II™ Fall Monitor

An intelligent choice in fall prevention, with Posey's Lifetime Warranty.

Lifetime Warranty on Alarm Unit

Plus one-year warranty on reusable sensors for bed and chair. Assures quality and long-lasting dependability.



Voice Recording

Records a personalized message in patient's own language, which reduces agitation, minimizes noise, and extends the caregiver's presence.



4 Alarm Modes

Multiple options include voice only (with nurse call notification), tone only, voice plus tone, or mute (with nurse call), allowing you to select the most appropriate type of alarm notification for individual patients/residents.

Hold Feature

Maintains setting for up to 6 hours while the patient/resident is away from bed or chair. The patient can be moved for therapy, social functions, meals, etc., without caregivers needing to reset the alarm. This ensures continuity of care and monitoring despite multiple caregivers.



10 Alarm Tones

Ten different tones – including a selection of familiar tunes – let you identify patients individually and differentiate from other alarms in the facility.

Light Indicators

Easy-to-see indicators for normal operation, hold, and low battery assure normal operation day or night, without having to disturb the patient/resident.

Easy Set-up and Operation

Sensor cables are pre-connected for quick set-up. Only one simple-to-use phone jack connection is required.

Variety of Sensors Assures Versatility

May be used with chair pad, bed pad, or chair belt sensors. Bed sensors include above-mattress or below-mattress pads to accommodate most institutional bed and mattress styles, including mattress overlays. Under-mattress sensor prevents tampering, minimizes false alarms. Available options include single-patient disposable sensor or reusable sensor, offering cost-effective solutions to both short and long-term stay situations.

Highly Visible LCD Display

Settings – alarm mode, volume, tone and variable delay – are continuously displayed for quick verification.

Optional Nurse Call Interface

The nurse-call interface can be used in conjunction with the in-room audible alarm, or the in-room alarm may be silenced to allow undisturbed sleep.

Adjustable Alarm Delay

Delay of up to 10 seconds (user-adjustable) before alarm sounds reduces false alarms and accommodates restless sleepers.

Adjustable Alarm Volume

Five volume settings allow the caregiver to set the alarm volume to suit facility and patient needs. The volume ranges from 80 to 103 decibels.

Battery-Operated

The unit runs on 4 "C" batteries to allow portability. Battery life is estimated at 30 days of normal use. Low battery warning light gives ample notice that a battery change is needed.



Multiple Mounting Options

Easily attaches to bed, chair, wall, or IV pole, accommodating a variety of uses.

Fail-Safe Sensor Alarm

If the sensor is removed from the unit while armed, the monitor will alarm to alert the caregiver of a possible monitoring interruption.

Sitter II
The fall monitor
that thinks for itself



SITTER II SPECIFICATIONS

Order Numbers	Sitter II - No Sensors	8281	Gel Foam Alarm Cushion	6509	Bed Bracket, standard	8276
	Chair Belt System	8280B	Sitter II Reusable Sensor	8282H/L	Bed Bracket, large	8277
	Sitter II complete, Chair Sensor	8280C	Above-Mattress Bed Sensor	8283	Sitter II Cradle	8279
	Gel Foam Cushion System	8280G	Sitter II Chair Pad	8355	IV Pole Mount	8288
	Below-Mattress System, heavy	8280H	Chair Belt Sensor	8360	Wheelchair Bracket	8289
	Below-Mattress System, light	8280L	Nurse Call Cable	8284	"Y" Adapter	8235S
	Above-Mattress System	8280M				
Size	Alarm Unit 4 1/2" wide x 7 1/2" long x 2" deep		Maintenance	Virtually maintenance-free		
Weight	1.3 lbs. with batteries; 0.7 lbs. without batteries		Warranty	Lifetime Warranty on Alarm Unit, parts and labor One-year Warranty on Chair Belt Sensor One-year Warranty on Below-Mattress Bed Pad Sensor 30-day Warranty on Above-Mattress Bed Pad Sensor		
Power Supply	4 alkaline, 1.5 volt "C" cell batteries					

Lifetime Warranty – The Posey Company warrants, to the original purchaser, this alarm unit to be free from defects in materials and workmanship for the life of the unit. Warranty void if unit is opened, damaged, or submerged.

Free 30-day trial. Call 1-800-447-6739

J.T. Posey Company, 5635 Peck Road, Arcadia, CA 91006-0020 • Phone 800-447-6739 • Fax 800-767-3933 • www.posey.com

©2000 J.T. Posey Company M1065 040400